



Additional options for RCN Shower and Commode Chair range SC 200 and SCC 200



BS Ergonomically designed monobloc seat with toilet and hygiene opening, made of solid plastic



PUS Plain soft seat made of PUR-foam



PS Soft seat made of PUR-foam with toilet opening



PSE Insert for soft seat PS



PP Soft seat with toilet and hygiene opening, made of PUR-foam



PPE Insert for soft seat PP



PZ Soft split seat made of PUR-foam



SP Seat cushion on plastic plate, with foam core and cover made of coated fabric
These cushions must be taken off prior to showering



SB Safety-bar




KA Drop-down armrest on the right, by request on the left



HV Height adjustment to adapt seat height and clearance to 5 levels (7 levels on chairs without ASA)



AP Armrest pad, made of PUR-foam




UR Splash protector made of PUR-foam, fits all PP seats



DRS Anti-slide, flips down and can be taken off



SG Safety-belt, continuous adjustment, easy to attach on frame



TK Second pail, handle upon request



NH Slant seat, 3 cm lower in back



AKA Anti-tip frame for very agitated users



FT Slide-out, anti-tip footrest



MFT Roll-along footrest



ER Ø100 mm single-wheel never rust caster - locking