

# The RCN product line for children



[www.rcn-medizin.com](http://www.rcn-medizin.com)



# Why RCN?


## RCN products for children

Especially impaired children and adolescents need special care and attention. This asks for healthcare equipment that makes the daily caring easier for all concerned.

RCN offers a broad range of products that are perfectly suited to respond to the needs of these younger people. RCN has been a manufacturer of quality healthcare equipment for over 30 years, striving to provide high quality products and to offer the market better solutions through innovation and by staying in touch with its customers.

With our products we focus on the field of personal hygiene which often requires a very special support - whether while showering, bathing, washing or visits to the toilet. As a result of the use of high-quality materials our products convince through their longevity and no-

rust guarantee. Products are made of RCN proprietary high impact plastic which is easy to clean and quickly warms up to body temperature. Moreover, we offer a great variety of equipment and color options to further adapt the product to individual needs and desires. Also RCN's proprietary modular system is ideally suited to respond to very specific product requirements.



We gladly realize individual requests - just contact us!



# Because expertise and quality count.

If one of our standard products does not fit your personal requirements we can easily make custom-built products that best respond to the individual needs of our little customers.

Discover the RCN product line for children on the following pages. If you have any questions or need further information, please let us know and do not hesitate to contact us.


## Our expertise - your gain

### RCN products:

- Are field tested by healthcare professionals
- Require very little maintenance
- Can be economically adapted to individual requirements
- Are tested by the TÜV lab and offer exceptional longevity

### Our promise to you:

- Short delivery times
- Specific product solutions
- Great customer service
- Lowest cost of ownership




# The Multi-Purpose Shower and Commode Chair



The KSC 180, the best solution for showering, toileting and transfer.

## KSC 180 FT RP SP

- Height adjustment to 7 levels, allows to increase the seat height while the child is growing
- Flat plastic seat with toilet opening
- Slide-out footrest FT
- Back support made of Textilene® fabric
- Removable backrest cushion RP and seat cushion SP with foam core and cover made of Textilene® fabric
- Ø 75 mm twin-wheel swivel casters with lock
- Pail with inside scale
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- supports up to 150kg
- 3 year warranty
- 5 year warranty on frame



Seat height adjustment from 37,5 to 52,5 cm

# Options for the KSC 180



Backrest cushion with foam core and cover made of Textilene® fabric RP, optional



Seat cushion on plastic plate, with foam core and cover made of Textilene® fabric SP, optional



Backrest support made of Textilene® fabric, standard



Flat plastic seat with toilet opening and pail with inside scale, standard



Slide-out footrest FT, optional



Swing-out armrest ASA, optional



Safety bar SB, optional

**New optional soft seat PPK, with toilet and hygiene opening; heat welded foam cushion for optimal comfort. Available in Grey.**

## Colors for seat cushion and backrest support made of Textilene® fabric "Sunsure"



white



navy



blue



suncastr  
blue



wine-  
berry



laguna



yellow



confetti

Model	Exterior dimensions			Seat			Lower leg length	Max seat clearance	Capacity	Weight
	Depth	Width	Height	Depth	Width	Height				
KSC 180 FT RP SP	55,5	46	84-99,5	41,5*	35	37,5-52,5*	25-40	29,5-44,5**	150	12,5

\*without cushion

\*\*Please note also the maximum width of 34cm underneath the seat

All information in cm resp. kg.

All information subject to change without prior notice!

# The Wheelchair for children DR 400 Mini



## DR 400 Mini SB

The wheelchair for children - suitable for self-drivers

- 22" wheels
- PUR-foam padded armrests
- Seat- and Backrest Cushion with foam core and cover made of water repellent coated fabric
- Drop-out footrest
- 2 wheel locks
- Safety bar **SB**, padded with PUR-foam
- Push-bar
- Clothes protector on both sides
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 120 kg
- 3 year warranty
- 5 year warranty on frame



## Leg support BSM

- Leg support **BSM** optional for DR 400 Mini, adjustable in 6 angle positions

The DR 400 Mini is not suitable for use in the shower.

### Coated fabric for seat- and backrest upholstery



saffron cherry sky pine charcoal

### PUR-foam for padded armrests and safety bar



graphite blue cobalt laguna

## Leg support BSM



Leg support adjustable in 6 angle positions.



Seat- and Backrest Cushion with foam core and cover made of water repellent coated fabric



PUR-foam padded armrests



Drop-out footrest



Clothes protector on both sides



2 wheel locks



Push-bar



Safety bar **SB**, padded with PUR-foam, optional



Leg support **BSM**, optional

Model	Exterior dimensions			Seat			Lower leg length	Capacity	Weight
	Depth	Width	Height	Depth	Width	Height			
DR 400 Mini SB	74	63,5	101	40	36	48	34	120	23

All information in cm resp. kg

All information subject to change without prior notice!

# The Shower-Wheelchair for children DR 100 Mini S



Backrest support made of Textilene® fabric "Sunsure"



white



navy



blue



suncastr  
blue



wine-  
berry



laguna



yellow



confetti

## DR 100 Mini S

The DR 100 Mini S is perfectly suitable for shower, pool or hydrotherapy. The high stability even on wet floors and the no-rust guarantee make this wheelchair an optimal tool for children and adolescents that want to move independently.

- 22" wheels
- 2 single-wheel swivel front casters, Ø 100 mm
- Plastic seat with drain holes
- Backrest made of Textilene® fabric
- Slide-out footrest FT
- Removable armrest on both sides 2 ASA
- Frame made of RCN proprietary high impact plastic
- Supports up to 120 kg
- 3 year warranty
- 5 year warranty on frame



Backrest made of Textilene® fabric



Plastic seat with drain holes



22" wheels



2 wheel locks



single-wheel swivel casters at the front, Ø 100 mm



Slide-out footrest FT



Removable armrest on both sides 2 ASA



Model	Exterior dimensions			Seat			Lower leg length	Capacity	Weight
	Depth	Width	Height	Depth	Width	Height			
DR 100 Mini S	78	64	76	38	36	46,5	29	120	18

All information in cm resp. kg

All information subject to change without prior notice!

# The Walker for children



Backrest support made of Textilene® fabric "Sunsure"



white



navy



blue



suncastr  
blue



wine-  
berry



laguna



yellow



confetti

## GW 120 Mini

The RCN Walker offers safety in movement like no other walking aid. Its equipment supports children's urge to move within secure conditions. Due to the easy-to-use height adjustment, this Walker is suitable for children

up to a height of 1.55 meters. The anti-tip, all around frame and the integrated seat together with the safety-belt greatly reduces fear of falling and protects against falling.



Anti-tip safety-frame, standard



Safety-bar, standard



Seat cushion with foam core and cover made of water repellent coated fabric, standard



Height adjustment to 7 positions, standard



Tote bag, standard



Backrest support made of Textilene® fabric, standard



6 Ø 75 mm twin-wheel swivel casters - 4 locking, standard



Sliding footrest GFT (can easily be taken off), optional



Drip stand GIH, optional



Safety-belt pad GGP, optional



Lap tray GSB, optional



Wall protectors made of PUR-foam AFS, optional

Model	Exterior dimensions			Seat			Grip height	Capacity	Weight
	Depth	Width	Height	Depth	Width	Height			
GW 120 Mini	89	72	88,5-103,5	26	38,5	39-54	59,5-74,5	150	16

All information in cm resp. kg

All information subject to change without prior notice!

# Bathing Support and Shower Gurneys



Continuously adjustable back- and legrest



## KDL 100

The series of KDL bathing and shower gurneys are ideally suitable for use in shower and bathtubs. They provide support for childrens who can not sit upright without assistance. The KDL bathing gurneys impress through their continuously adjustable back- and legrest and the very sturdy and comfortable sling deck made of Textilene® fabric. All materials are easy to clean and quickly adapt to body temperature, thus supporting a high-quality and pleasant care.

- Frame made of RCN proprietary high impact plastic
- Continuously adjustable back- and leg rest
- Robust sling deck made of Textilene® fabric
- Guaranteed NO RUST
- UV-resistant and compatible with most disinfectants
- 3 year warranty
- 5 year warranty on frame

Cover made of Textilene® fabric "Sunsure"



white



navy



blue



suncastr  
blue



wine-  
berry



laguna



yellow



confetti

Model	Exterior dimensions			Length of reclining areas			Capacity	Weight
	Length	Width	Height	Backrest	Seat	Legrest		
KDL 100	110	48,5	27	40	40	30	75	7
KDL 200	140	48,5	27	50	50	40	100	9

All information in cm resp. kg

All information subject to change without prior notice!

## UDL 500

The RCN shower gurneys are ideal for transporting, showering and bathing children and adolescents. The frame made of RCN proprietary high impact plastic guarantees long durability and a low tare weight with simultaneous high stability. The padded soft tub-liner offers comfort, is easy to clean and is compatible with most disinfectants. Closing of the drain opening allows to keep warm water in the tub, offering a bathtub feel.

- Padded soft tub-liner
- Drop-down side rails
- Large drain hose
- UV-resistant and compatible with most disinfectants
- Guaranteed NO RUST
- 4 single-wheel swivel casters (plastic/stainless steel), all locking
- Supports up to 150 kg
- For people up to 146 cm in height
- Shoulder width up to 48 cm
- 3 year warranty
- 5 year warranty on frame made of RCN proprietary high impact plastic




Drain hose, 90 cm



Closable drain opening



6 positions backrest adjustment **VR**, optional



Wedge pillow **DK**, optional



Drop-down side rails

Model	Exterior dimensions			Reclining area			Capacity	Weight
	Length	Width	Height	Up to body length	Up to shoulder width	Height		
UDL 500	164	75	96	146	48	80	150	41
UDL 750	180	80	96	162	53	80	150	46

All information in cm resp. kg

All information subject to change without prior notice!

# Height adjustable Shower Gurneys

## The UDW series

The UDW Shower Gurneys from RCN are ideal for transporting and showering young people in need of care. These gurneys are made of a powder coated stainless steel frame and are equipped with a padded soft tub-liner. More features, such as drop-down side and foot-end rails, fold down head-end rail, height adjustment,

Anti- / Trendelenburg positioning and an optional continuously adjustable backrest, make this high-quality gurney a very special nursing aid. The drain plug in connection with a 90 cm drain hose ensures that water drains off in a controlled manner.

- Padded soft tub-liner
- Drop-down side- and foot-end rails
- Head-end rail can be put in a lateral position to extend the length of the deck by 15 cm
- Anti- / Trendelenburg positioning via manually activated pneumatic-spring
- Height adjustable from 49 to 89 cm via footpedals on both sides
- Closable drain opening and drain hose, 90 cm
- Frame made of powder-coated stainless steel
- Supports up to 150 kg
- For people up to 162 / 177 cm in height
- Shoulder width up to 54 cm
- 8 twin-wheel swivel casters Ø 100 mm (4 locking) offer easy transport and great mobility
- Guaranteed NO RUST
- Easy to disinfect and to clean
- 3 year warranty
- 5 year warranty on frame



Picture: UDW 650

Drop-down side- and foot-end rails and and head-end rail extension



Height adjustment from 49 to 89 cm

Anti- / Trendelenburg positioning



Guaranteed NO RUST!



Continuously adjustable backrest **VRW**, optional

Picture: **UDW 650 VRW**

Model	Exterior dimensions			Reclining area			Capacity	Weight
	Length	Width	Height	Up to body length	Up to shoulder width	Height		
UDW 650	177	71	70-110	162/177*	54	49-89	150	69
UDW 650 VRW	177	71	70-110	162/177*	54	49-89	150	69

All information in cm resp. kg

\*With extended head rail

All information subject to change without prior notice!

# The baby changing unit WK 07



The changing unit WK 07 includes the bathtub cart **BWR 200** as a standard feature.

## WK 07

The mobile changing unit WK 07 makes the care of babies much more easier. The padded changing table has a comfortable working height of 97 cm which can be customized upon request. The lateral integrated cabinet provides enough space for all care utensils

needed and in combination with the mobile bathtub cart you have everything within reach. All materials used are compatible with most disinfectants and resistant to urine and blood.



# Only the best for the little ones

## WK 07

- Frame made of RCN proprietary high impact plastic
- Padded changing surface and three side padded rail with 30mm foam core and cover made of coated fabric
- Integrated cabinet with 1 storage compartment and 2 drawers made of scratch resistant RESOPAL® HPL plates
- 6 twin-wheel swivel casters, Ø 75 mm, 4 locking
- Guaranteed NO RUST
- 3 year warranty
- 5 year warranty on frame

with **BWR 200**:

- Mobile bathtub cart
- Frame made of RCN proprietary high impact plastic
- Bathtub color blue, including 70cm drain hose
- 4 twin-wheel swivel casters, Ø 75 mm, all locking



The mobile bathtub stores easily underneath the changing deck, saving space.




The Bathtub Cart BWR 200 is available independently.

### Upholstery changing unit:



saffron cherry sky pine charcoal

### Cabinet made of RESOPAL® HPL:



leafgreen gentian blue signal red clementine orange broom yellow alabaster white white

Model: WK 07	Width	Depth	Height	Weight
Exterior dimensions	101	86	117	74
Surface of changing pad	82,5	76	97*	-
Inner dimension drawer	24	58	16	-
BWR 200	76	45,5	87,5	8,6

All information in cm resp. kg. All information subject to change without prior notice!

\*Height of work surface

# The changing unit WK 07 with heat lamp

## The RCN changing unit WK 07 WL

The targeted and controlled heat supply IW 980 AGU has been especially developed for a targeted heat supply. It protects babies against heat losses and helps to maintain their body temperature.

In connection with a skin temperature control function the heat lamp ensures a constant baby's skin temperature by automatically adjusting the heat emission to physiological and environmental conditions. This is achieved due to a microprocessor which checks the skin temperature 10 times per second. Thus babies are wonderfully relaxed while being taken care of.



Heat lamp easily swings away side-ways.



## WK 07 WL

All features just like WK 07  
but with additional:

- Certified heat lamp holder
- Heat lamp unit IW 980 AGU
- Skin temperature control
- DuoSense skin temperature sensor
- SunSpot sensor protection
- IV holder

# Children want to enjoy nature



## The Beach Mobil SM 700 FBS

The Beach Mobil SM 700 FBS made of white RCN proprietary high impact plastic - break-proof, guaranteed NO RUST, UV-resistant. Thanks to its low pressure balloon tires the chair rolls easily over the most soft, sandy beaches and other uneven surfaces.

Seat- and backrest cushion made with foam core and Textilene® fabric are comfortable for hours of beach enjoyment at the seaside where standard wheel chairs cannot go.

The Beach Mobil can be knocked down into three sections – no tools needed – making transport to and from the beach more amenable.

This Beach Mobil lets wheel chair bound young people enjoy the pleasures of a day on the beach!

- Hand operated brake/lock
- Tote bag
- NO RUST
- Salt water resistant
- Ideal for sandy beaches and other uneven and soft surfaces



Options:

- Safety bar **SBS**
- Parasol with holder **SO**

Seat- and backrest cushion made of Textilene® fabric "Sunsure"



white



navy



blue



suncastr  
blue



wine-  
berry



laguna



yellow






confetti

SM 700 FBS	Length	Width	Height	Weight
Overall dimensions	110	104	111	35

All information in cm resp. kg

All information subject to change without prior notice!



**RCN - the best address for  
Healthcare Equipment, Quality and Service!**

**RCN Medizin- und Rehattechnik GmbH**

Hauptstraße 4 · D-55471 Sargenroth

Tel. +49 (0) 6761-9197-0 · Fax +49 (0) 6761-9197-70

info@rcn-medizin.de

[www.rcn-medizin.com](http://www.rcn-medizin.com)

01/2018

