

Datasheet: AT 2 MRT

MRI compatible Step Stool

Price ex works Sargenroth:

Net list price [EUR]: 327,00

Technical data [mm/kg]:

Height: 380

Width: 400

Depth: 660

Weight: 7,5

Height of the first step: 200

Weight capacity: 150

Subject to technical changes!



Date: 01.01.2024

Features:

- Fulfills all safety requirements for stations up to 3 Tesla - MR CONDITIONAL
- Step stool made of solid plastic, double-stage
- No-slip strips on the tread surface
- Special anti-slip legs
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Easy to clean
- Supports up to 150 kg
- 3 year warranty
- 5 year warranty on frame
- For stations higher than 3 Tesla we offer customer specific product solutions, even all plastic equipment if required