

## Price ex works Sargenroth:

Net list price [EUR]: 659,00

## Technical data [mm/kg]:

Height:	875
Width:	455
Length:	775
Weight:	8,6
Length of bathtub (inside):	580
Width of bathtub (inside):	260
Height of bathtub (inside):	200
Subject to technical changes!	



Date: 01.01.2024

## Features:

- Different height by request without extra charge
- Bathtub Cart - frame made of RCN proprietary high impact plastic
- Bathtub color blue, including drain hose (hose length 840 mm)
- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm, all locking
- UV-resistant and compatible with most disinfectants
- Easy to clean
- 3 year warranty
- 5 year warranty on frame

## Color selection bathtub



light blue    gray