

Datasheet: DES 200

## High chair for X-ray defecography examination

### Price ex works Sargenroth:

Net list price [EUR]: 2.097,00

### Technical data [mm/kg]:

Height:	1005
Width:	450
Depth:	1210
Weight:	22
Seat height:	800
Height of the first step:	200
Height of the second step:	400
Weight capacity:	150

Subject to technical changes!



Date: 01.01.2024

### Features:

- For x-ray examination and diagnosis of anorectal dysfunction
- No materials blocking the X-rays ensures high precision pictures of the affected area
- Backrest, padded with PUR-foam
- Oval pail and lid included
- Fitting step stool (single or double level)
- Seat height and step height according to individual specifications
- Steps can be easily removed from the chair without tools
- Special anti-slip legs
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Easy to clean
- 3 year warranty
- 5 year warranty on frame

### Color selection PUR-foam

