

# Datasheet: DH 100 Shower Chair

## Price ex works Sargenroth:

Net list price [EUR]: 439,00

## Technical data [mm/kg]:

Height:	815
Width:	545
Depth:	545
Weight:	8,7
Seat height:	530
Seat width:	440
Seat depth:	500
Clearance under seat:	450
Weight capacity:	150
Subject to technical changes!	



Date: 01.01.2024

## Features:

- Armrests
- Backrest
- Ergonomically designed toilet seat, lifts up and can be removed
- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm, all locking
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Easy to clean
- 3 year warranty
- 5 year warranty on frame
- Guaranteed NO RUST