

# Datasheet: DH 100 PP

## Shower Chair

### Price ex works Sargenroth:

Net list price [EUR]: 642,00

### Technical data [mm/kg]:

|                       |      |
|-----------------------|------|
| Height:               | 815  |
| Width:                | 545  |
| Depth:                | 545  |
| Weight:               | 10,1 |
| Seat height:          | 550  |
| Seat width:           | 440  |
| Seat depth:           | 500  |
| Clearance under seat: | 450  |
| Weight capacity:      | 150  |

Subject to technical changes!



Date: 01.01.2026

### Features:

- Armrests
- Backrest
- Soft seat with toilet and hygiene opening made of PUR-foam
- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm, all locking
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Easy to clean
- 3 year warranty
- 5 year warranty on frame
- Guaranteed NO RUST
- Option for more comfort: insert for soft seat (PPE)

### Color selection PUR-foam



graphite blue cobalt