

Datasheet: DH 49

Shower Stool

Price ex works Sargenroth:

Net list price [EUR]: 119,00

Technical data [mm/kg]:

Height:	480
Width:	420
Depth:	420
Weight:	3,4
Seat height:	470
Seat width:	420
Seat depth:	420
Weight capacity:	200

Subject to technical changes!



Date: 01.01.2024

Features:

- Ergonomically designed seat with central hygiene opening
- Special anti-slip legs
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 200 kg
- Easy to clean
- 3 year warranty
- 5 year warranty on frame
- Guaranteed NO RUST
- Option for more comfort: seat pad made of PUR-foam (PA)