

# Datasheet: DH 49 DLR Shower Stool

## Price ex works Sargenroth:

Net list price [EUR]: 223,00

## Technical data [mm/kg]:

|                  |     |
|------------------|-----|
| Height:          | 485 |
| Width:           | 520 |
| Depth:           | 520 |
| Weight:          | 5,2 |
| Seat height:     | 470 |
| Seat width:      | 420 |
| Seat depth:      | 420 |
| Weight capacity: | 150 |

Subject to technical changes!



Date: 01.01.2024

## Features:

- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm, all locking
- Ergonomically designed seat with central hygiene opening
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Easy to clean
- 3 year warranty
- 5 year warranty on frame
- Guaranteed NO RUST
- Option for more comfort: seat pad made of PUR-foam (PA)