

Datasheet: DH 49 DLR PA MRT

MRI compatible stool

Price ex works Sargenroth:

Net list price [EUR]: 441,00

Technical data [mm/kg]:

Height:	500
Width:	520
Depth:	520
Weight:	6
Seat height:	485
Seat width:	420
Seat depth:	420
Weight capacity:	150

Subject to technical changes!



Date: 01.01.2024

Features:

- Fulfills all safety requirements for stations up to 3 Tesla - MR CONDITIONAL
- 4 twin-wheel swivel casters (plastic/brass), Ø 75 mm, all locking, non-magnetic
- Ergonomically designed seat with central hygiene opening
- Seat pad made of PUR-foam
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Easy to clean
- 3 year warranty
- 5 year warranty on frame
- For stations higher than 3 Tesla we offer customer specific product solutions, even all plastic equipment if required

Color selection PUR-foam

