

# Datasheet: DH 49 PA MRT

## MRI compatible stool

### Price ex works Sargenroth:

Net list price [EUR]: 213,00

### Technical data [mm/kg]:

Height:	495
Width:	420
Depth:	420
Weight:	4,7
Seat height:	485
Seat width:	420
Seat depth:	420
Weight capacity:	200

Subject to technical changes!



Date: 01.01.2024

### Features:

- Fulfills all safety requirements for stations up to 3 Tesla - MR CONDITIONAL
- Ergonomically designed seat with central hygiene opening
- Seat pad made of PUR-foam
- Special anti-slip legs
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 200 kg
- Easy to clean
- 3 year warranty
- 5 year warranty on frame
- For stations higher than 3 Tesla we offer customer specific product solutions, even all plastic equipment if required

### Color selection PUR-foam



graphite blue cobalt