

Datasheet: DH 49 PRS Shower Stool

Price ex works Sargenroth:

Net list price [EUR]: 268,00

Technical data [mm/kg]:

Height:	790
Width:	430
Depth:	550
Weight:	6,9
Seat height:	470
Seat width:	420
Seat depth:	420
Weight capacity:	150

Subject to technical changes!



Date: 01.01.2024

Features:

- Soft backrest made of PUR-foam, only in blue
- Ergonomically designed seat with central hygiene opening
- Special anti-slip legs
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Easy to clean
- 3 year warranty
- 5 year warranty on frame
- Guaranteed NO RUST
- Option for more comfort: seat pad made of PUR-foam (PA)