

# Datasheet: DH 49 RL

## Shower Stool

### Price ex works Sargenroth:

Net list price [EUR]: 181,00

### Technical data [mm/kg]:

Height:	660
Width:	430
Depth:	500
Weight:	4,9
Seat height:	470
Seat width:	420
Seat depth:	420
Weight capacity:	150

Subject to technical changes!



Date: 01.01.2024

### Features:

- Backrest, padded with PUR-foam
- Ergonomically designed seat with central hygiene opening
- Special anti-slip legs
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Easy to clean
- 3 year warranty
- 5 year warranty on frame
- Guaranteed NO RUST
- Option for more comfort: seat pad made of PUR-foam (PA)

### Color selection PUR-foam



graphite blue cobalt