

# Datasheet: DH 52

## Shower Stool

### Price ex works Sargenroth:

Net list price [EUR]: 265,00

### Technical data [mm/kg]:

|                  |     |
|------------------|-----|
| Height:          | 500 |
| Width:           | 500 |
| Depth:           | 340 |
| Weight:          | 7   |
| Seat height:     | 500 |
| Seat width:      | 500 |
| Seat depth:      | 330 |
| Weight capacity: | 150 |

Subject to technical changes!



Date: 01.01.2024

### Features:

- Plastic seat with handle on both sides
- Special anti-slip legs
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Easy to clean
- 3 year warranty
- 5 year warranty on frame
- Guaranteed NO RUST