

Datasheet: DH 80

Shower Stool

Price ex works Sargenroth:

Net list price [EUR]: 310,00

Technical data [mm/kg]:

Height:	685
Width:	540
Depth:	510
Weight:	7,1
Seat height:	510
Seat width:	440
Seat depth:	510
Clearance under seat:	430
Weight capacity:	150
Subject to technical changes!	



Date: 01.01.2024

Features:

- Armrests
- Ergonomically designed toilet seat, lifts up and can be removed
- Special anti-slip legs
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Easy to clean
- 3 year warranty
- 5 year warranty on frame
- Guaranteed NO RUST