

# Datasheet: DH 80 PP

## Shower Stool

### Price ex works Sargenroth:

Net list price [EUR]: 459,00

### Technical data [mm/kg]:

Height:	685
Width:	540
Depth:	520
Weight:	8,3
Seat height:	530
Seat width:	440
Seat depth:	520
Clearance under seat:	430
Weight capacity:	150
Subject to technical changes!	



Date: 01.01.2024

### Features:

- Armrests
- Soft seat with toilet and hygiene opening made of PUR-foam
- Special anti-slip legs
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Easy to clean
- 3 year warranty
- 5 year warranty on frame
- Guaranteed NO RUST
- Option for more comfort: insert for soft seat (PPE)

### Color selection PUR-foam



graphite blue cobalt