

Datasheet: DH 90

Shower Stool

Price ex works Sargenroth:

Net list price [EUR]: 381,00

Technical data [mm/kg]:

Height:	795
Width:	540
Depth:	545
Weight:	7,6
Seat height:	510
Seat width:	440
Seat depth:	500
Clearance under seat:	430
Weight capacity:	150

Subject to technical changes!



Date: 01.01.2024

Features:

- Backrest
- Armrests
- Ergonomically designed toilet seat, lifts up and can be removed
- Special anti-slip legs
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Easy to clean
- 3 year warranty
- 5 year warranty on frame
- Guaranteed NO RUST