

Datasheet: DH 90 PP

Shower Stool

Price ex works Sargenroth:

Net list price [EUR]: 530,00

Technical data [mm/kg]:

Height:	795
Width:	540
Depth:	545
Weight:	9,3
Seat height:	530
Seat width:	440
Seat depth:	500
Clearance under seat:	430
Weight capacity:	150
Subject to technical changes!	



Date: 01.01.2024

Features:

- Backrest
- Armrests
- Soft seat with toilet and hygiene opening made of PUR-foam
- Special anti-slip legs
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Easy to clean
- 3 year warranty
- 5 year warranty on frame
- Guaranteed NO RUST

Color selection PUR-foam



graphite blue cobalt