

Datasheet: DR 100 K

All-plastic Transfer Chair

Price ex works Sargenroth:

Net list price [EUR]: 3.363,00

Technical data [mm/kg]:

Height:	1075
Width:	705
Depth:	855
Weight:	27
Seat height:	520
Seat width:	475
Seat depth:	440
Weight capacity:	120
Lower leg length:	460
Subject to technical changes!	



Date: 01.01.2024

Features:

- Completely made of plastic, without any metal parts
- Suitable for all MRI systems and electric security gates (does not activate an alarm) - MR SAFE
- Suitable for self-drivers
- Seat- and backrest pad made of PUR-foam, only available in black
- Clothes protector on both sides
- Angled armrests, PUR-foam padded
- Drop-out footrest
- Push-bar
- Puncture-proof tires
- 24" wheels
- 2 single-wheel swivel casters (all-plastic) at the front, Ø 125 mm
- 2 wheel locks
- Supports up to 120 kg
- UV-resistant and compatible with most disinfectants
- 2 year warranty