

Datasheet: DR 100 MRT

MRI compatible Transport Chair up to 3 Tesla

Price ex works Sargenroth:

Net list price [EUR]: 2.618,00

Technical data [mm/kg]:

Height:	1075
Width:	725
Depth:	800
Weight:	32
Seat height:	570
Seat width:	440
Seat depth:	450
Lower leg length:	400
Width without handrail:	675
Weight capacity:	150

Subject to technical changes!



Date: 01.01.2024

Features:

- Fulfills all safety requirements for stations up to 3 Tesla - MR CONDITIONAL
- Suitable for self-drivers
- 24" wheels with brass axles and plastic bearings
- Puncture-proof tires
- Swing-out armrest on both sides
- Slide-out footrest
- Soft backrest made of PUR-foam, removable for easy cleaning and access
- Soft seat made of PUR-foam
- Push-bar
- Angled armrests
- With slightly inclined seat (3 degrees)
- 2 wheel locks
- 2 single-wheel swivel casters (all-plastic) at the front, Ø 125 mm
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- 3 year warranty
- 5 year warranty on frame
- For stations higher than 3 Tesla we offer customer specific product solutions, even all plastic equipment if required

Color selection PUR-foam



graphite blue cobalt