

# Datasheet: DR 100 Mini S Shower-Wheelchair

## Price ex works Sargenroth:

Net list price [EUR]: 1.993,00

## Technical data [mm/kg]:

Height:	760
Width:	640
Depth:	750
Weight:	18
Seat height:	470
Seat width:	360
Seat depth:	380
Weight capacity:	120
Lower leg length:	350
Subject to technical changes!	



Date: 01.01.2024

## Features:

- Suitable for self-drivers
- Plastic seat with drain holes
- Backrest made of "Sunsure" fabric
- Slide-out footrest
- Removable armrest on both sides
- Puncture-proof tires
- 22" wheels
- 2 wheel locks
- 2 twin-wheel swivel casters (plastic/stainless steel) at the front, Ø 75 mm, non-locking
- Easy to clean
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- 3 year warranty
- 5 year warranty on frame

## Color selection Sunsure fabric

