

# Datasheet: DR 100 OS PPG XXL Shower-Wheelchair

## Price ex works Sargenroth:

Net list price [EUR]: 3.062,00

## Technical data [mm/kg]:

Height:	1010
Width:	840-880
Depth:	790
Weight:	37
Seat height:	540
Seat width:	560-600
Seat depth:	450
Weight capacity:	200
Lower leg length:	375

Subject to technical changes!



Date: 01.01.2024

## Features:

- 24" wheels
- Swing-out armrest on both sides
- Slide-out, roll-along footrest
- Soft seat with toilet and hygiene opening made of PUR-foam, only available in blue
- Soft backrest made of PUR-foam, removable for easy cleaning and access
- 2 wheel locks
- 2 single-wheel swivel casters (plastic/stainless steel) at the front, Ø 100 mm
- Footrest with 4 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm
- Suitable for self-drivers
- Does NOT roll over toilet
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 200 kg
- Guaranteed NO RUST
- 3 year warranty
- 5 year warranty on frame