

Datasheet: DR 200 S

Shower- / Commode Wheelchair

Price ex works Sargenroth:

Net list price [EUR]: 1.736,00

Technical data [mm/kg]:

Height:	1070
Width:	720
Depth:	790
Weight:	22
Seat height:	550
Seat width:	440
Seat depth:	450
Clearance under seat:	425
Lower leg length:	370
Weight capacity:	150

Subject to technical changes!



Date: 01.01.2024

Features:

- 24" wheels
- Puncture-proof tires
- 2 single-wheel swivel casters (plastic/stainless steel) at the front, Ø 125 mm
- Swing-out armrest on both sides
- PUR-foam padded armrests
- Slide-out footrest
- Monobloc seat with toilet and hygiene opening, made of solid plastic
- Monobloc backrest with push handles, made of solid plastic
- 2 wheel locks
- Pail with inside scale
- Rolls over toilet
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Guaranteed NO RUST
- Easy to clean
- 3 year warranty
- 5 year warranty on frame

Color selection plastic backrest / plastic seat



cobalt graphite red