

Datasheet: DR 250 OS PPG

XXL Shower- / Commode Wheelchair

Price ex works Sargenroth:

Net list price [EUR]: 2.695,00

Technical data [mm/kg]:

Height: 1010
Width: 840-880
Depth: 790
Weight: 39
Seat height: 540
Seat width: 560-600
Seat depth: 450
Weight capacity: 200
Lower leg length: 375
Subject to technical changes!



Date: 02.01.2021

Features:

- 24" wheels
- Swing-out armrest on both sides
- Slide-out, roll-along footrest
- Soft seat with toilet and hygiene opening made of PUR-foam, only available in blue
- Soft backrest made of PUR-foam, removable for easy cleaning and access
- 2 wheel locks
- 2 single-wheel swivel casters (plastic/stainless steel) at the front, Ø 100 mm
- Footrest with 4 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm
- Pail with inside scale
- Lid storage on backrest
- Suitable for self-drivers
- Does NOT roll over toilet
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 200 kg
- Guaranteed NO RUST
- 3 year warranty
- 5 year warranty on frame