

# Datasheet: DR 250 PPG

## Shower- / Commode Wheelchair

### Price ex works Sargenroth:

Net list price [EUR]: 2.162,00

### Technical data [mm/kg]:

Height:	1010
Width:	725
Depth:	800
Weight:	29
Seat height:	565
Seat width:	440
Seat depth:	450
Clearance under seat:	430
Lower leg length:	400
Width without handrail:	675
Weight capacity:	150

Subject to technical changes!



Date: 27.02.2024

### Features:

- 24" wheels
- Puncture-proof tires
- Swing-out armrest on both sides
- Slide-out footrest
- Soft backrest made of PUR-foam, removable for easy cleaning and access
- Soft seat with toilet and hygiene opening made of PUR-foam
- With slightly inclined seat (3 degrees)
- Suitable for self-drivers
- Rolls over toilet
- Pail with inside scale
- 2 wheel locks
- 2 single-wheel swivel casters (plastic/stainless steel) at the front, Ø 125 mm
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Guaranteed NO RUST
- 3 year warranty
- 5 year warranty on frame

### Color selection PUR-foam

