

Datasheet: DR 400

Transport Wheelchair

Price ex works Sargenroth:

Net list price [EUR]: 1.998,00

Technical data [mm/kg]:

Height:	1055
Width:	725
Depth:	800
Weight:	29
Seat height:	555
Seat width:	440
Seat depth:	400
Weight capacity:	150
Lower leg length:	350

Subject to technical changes!



Date: 01.01.2024

Features:

- 24" wheels
- Swing-out armrest on both sides
- Backrest cushion and seat cushion with foam core and cover made of coated fabric
- Slide-out footrest
- 2 wheel locks
- Push-bar
- Angled armrests, PUR-foam padded
- Clothes protector on both sides
- Suitable for self-drivers
- Puncture-proof tires
- 2 single-wheel swivel casters (plastic/stainless steel) at the front, Ø 125 mm
- Supports up to 150 kg
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- 3 year warranty
- 5 year warranty on frame

Color selection coated fabric

