

Datasheet: DR 400 Mini

Transport Wheelchair for children

Price ex works Sargenroth:

Net list price [EUR]: 1.217,00

Technical data [mm/kg]:

Height:	1010
Width:	635
Depth:	740
Weight:	23
Seat height:	480
Seat width:	360
Seat depth:	400
Weight capacity:	120
Lower leg length:	340

Subject to technical changes!



Date: 01.01.2024

Features:

- 22" wheels
- PUR-foam padded armrests
- Backrest cushion and seat cushion with foam core and cover made of coated fabric
- Drop-out footrest
- 2 wheel locks
- Push-bar
- Clothes protector on both sides
- Suitable for self-drivers
- Puncture-proof tires
- 2 single-wheel swivel casters (plastic/stainless steel) at the front, Ø 100 mm
- Supports up to 120 kg
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- 3 year warranty
- 5 year warranty on frame

Color selection coated fabric



Color selection PUR-foam

