

# Datasheet: DT 100 HM SP

## Shower Stool

### Price ex works Sargenroth:

Net list price [EUR]: 570,00

### Technical data [mm/kg]:

Height: 675-825

Width: 545

Depth: 510

Weight: 8,8

Seat height: 500-650

Seat width: 440

Seat depth: 500

Clearance under seat: 420-570

Weight capacity: 150

Subject to technical changes!



Date: 01.01.2024

### Features:

- All-in-One solution: Toilet Raiser, Rise-Up Aid, Shower Stool and Commode
- Armrests
- Ergonomically designed toilet seat, lifts up and can be removed
- With additional pail (for use as commode)
- With splash guard (for use as toilet raiser - simply put on top of toilet)
- Height adjustment to 7 different levels (25 mm steps)
- Seat cushion with foam core and cover made of coated fabric
- Special anti-slip legs
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Easy to clean
- 3 year warranty
- 5 year warranty on frame
- Guaranteed NO RUST

### Color selection coated fabric

