

Datasheet: DT 100 HV Shower Stool

Price ex works Sargenroth:

Net list price [EUR]: 417,00

Technical data [mm/kg]:

Height: 675-825

Width: 545

Depth: 510

Weight: 7,6

Seat height: 500-650

Seat width: 440

Seat depth: 500

Clearance under seat: 420-570

Weight capacity: 150

Subject to technical changes!



Date: 01.01.2024

Features:

- Toilet Raiser, Rise-Up Aid and Shower Stool, all in one
- Armrests
- Ergonomically designed toilet seat, lifts up and can be removed
- With splash guard (for use as toilet raiser - simply put on top of toilet)
- Height adjustment to 7 different levels (25 mm steps)
- Special anti-slip legs
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Easy to clean
- 3 year warranty
- 5 year warranty on frame
- Guaranteed NO RUST