

Datasheet: DWH 86

Shower Corner Stool

Price ex works Sargenroth:

Net list price [EUR]: 354,00

Technical data [mm/kg]:

Height:	615
Width:	570
Depth:	480
Weight:	6,7
Seat height:	510
Seat width:	435
Seat depth:	345
Weight capacity:	120

Subject to technical changes!



Date: 01.01.2024

Features:

- Large all plastic seat with open front
- Armrests
- Suction cups on front legs
- A smooth bathtub floor is necessary to ensure a firm and safe stand of the front legs with suction cups
- Only needs to be placed against two 90° walls
- Fits easily in even small shower stalls
- Best use of available space
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 120 kg
- Easy to clean
- Guaranteed NO RUST
- 3 year warranty
- 5 year warranty on frame
- Also available with shorter rear legs for a better positioning on the edge of the shower tray (from 10 cm edge height)
- Additional spacers are recommended from a distance of 23 cm (from shower corner to rear legs)
- Optional spacers: 55,00 EUR (set of 2 spacers)