

# Datasheet: GSS 800 S Rocking Chair

## Price ex works Sargenroth:

Net list price [EUR]: 1.626,00

## Technical data [mm/kg]:

Height:	970
Width:	665
Depth:	900
Weight:	31
Seat height:	500
Seat width:	540
Seat depth:	450
Depth with fold out footrest:	1110
Weight capacity:	120
Subject to technical changes!	



Date: 01.01.2024

## Features:

- Special suspension ensures rhythmic, sliding and nearly horizontal rocking movements
- Removable seat cushion (attached to seat plate using Velcro strips)
- Removable back cushion (fastened using snap fasteners)
- Backrest cushion and seat cushion with foam core and cover made of coated fabric
- Special non-tilting and stable legs
- Safety-bar
- PUR-foam padded armrests
- Footrest (folds upwards)
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Easy to clean
- 3 year warranty
- 5 year warranty on frame

## Color selection coated fabric



## Color selection PUR-foam

