

# Datasheet: GW 120 Maxi

## Maxi Walker

### Price ex works Sargenroth:

Net list price [EUR]: 1.459,00

### Technical data [mm/kg]:

Height: 1040-1190  
Width: 835  
Depth: 1125  
Weight: 30,5  
Seat height: 450-600  
Seat width: 500  
Seat depth: 340  
Grip height: 750-900  
Weight capacity: 200

Subject to technical changes!



Date: 01.01.2022

### Features:

- Seat cushion with foam core and cover made of coated fabric
- Backrest made of Textilene® fabric
- UV-resistant and compatible with most disinfectants
- Safety-bar and safety-belt
- Frame made of RCN proprietary high impact plastic
- 8 twin-wheel swivel casters, Ø 75 mm, 4 locking
- Height adjustment to 7 levels
- Tote bag made of Textilene® fabric
- Supports up to 200 kg
- Easy to clean
- 3 year warranty
- 5 year warranty on frame

### Color selection Textilene® „Sunsure“ fabric



### Color selection coated fabric

