

Datasheet: **HGW 100 BV**

Hygiene and Isolation Cart

Price ex works Sargenroth:

Net list price [EUR]: 1.810,00

Technical data [mm/kg]:

Height:	1385
Width:	975
Depth:	570
Weight:	38
Height of top plate:	865
Height of drawer:	165
Width of drawer:	530
Depth of drawer:	360

Subject to technical changes!



Date: 01.01.2026

Features:

- 10 small flip-out containers
- All plastic top plate
- 3 transparent drawers, each with a load capacity of 8 kg
- Liquid dispenser on slide rail (ESP), standard for half liter bottles (for liter bottles upon request)
- Drip tray for liquid dispenser
- Glove dispenser on slide rail (HSP)
- Holder for sharps disposable container
- Side-mounted waste container, 35 liter
- Slide-out writing platform
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm, 2 locking
- Easy to clean
- 3 year warranty
- 5 year warranty on frame