

Datasheet: **HGW 100 F**

Hygiene and Isolation Cart

Price ex works Sargenroth:

Net list price [EUR]: 769,00

Technical data [mm/kg]:

Height:	1230
Width:	700
Depth:	410
Weight:	19
Height of top plate:	865
Height of drawer:	165
Width of drawer:	530
Depth of drawer:	360

Subject to technical changes!



Date: 01.01.2026

Features:

- Holder for 0.5 l bottles with dosing pump on sliding rail
- Glove dispenser on slide rail (HSP)
- All plastic top plate
- 3 transparent drawers, each with a load capacity of 8 kg
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm, 2 locking
- Easy to clean
- 3 year warranty
- 5 year warranty on frame
- Delivery without pictured contents