

Datasheet: HGW 100 S

Hygiene and Isolation Cart

Price ex works Sargenroth:

Net list price [EUR]: 954,00

Technical data [mm/kg]:

Height: 1340

Width: 975

Depth: 410

Weight: 26

Height of top plate: 865

Height of drawer: 165

Width of drawer: 530

Depth of drawer: 360

Subject to technical changes!



Date: 01.01.2024

Features:

- Side-mounted waste container, 35 liter
- Liquid dispenser on slide rail (ESP), standard for half liter bottles (for liter bottles upon request)
- Triple glove dispenser on slide rail
- All plastic top plate
- 3 transparent drawers, each with a load capacity of 8 kg
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm, 2 locking
- Easy to clean
- 3 year warranty
- 5 year warranty on frame
- Delivery without pictured contents