

Datasheet: IH PH HV MRT

MRI compatible drip stand

Price ex works Sargenroth:

Net list price [EUR]: 508,00

Technical data [mm/kg]:

Height: 1530-1905

Width: 660

Depth: 660

Height of pump holder (in center): 1100

Diameter of pump holder: 32

Weight: 7

Weight capacity per pump holder: 9

Weight capacity per hook (cross above): 2

Total weight capacity of the cross: 4

Subject to technical changes!



Date: 01.01.2024

Features:

- Fulfills all safety requirements for stations up to 3 Tesla - MR CONDITIONAL
- Pump holder
- Height adjustment to 6 different levels (75 mm steps)
- 4 twin-wheel swivel casters (plastic/brass), Ø 75 mm, 2 locking, non-magnetic
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- 3 year warranty
- 5 year warranty on frame
- For stations higher than 3 Tesla we offer customer specific product solutions, even all plastic equipment if required