

Datasheet: KDL 200

Bathing support for children

Price ex works Sargenroth:

Net list price [EUR]: 1.234,00

Technical data [mm/kg]:

Height:	270
Width:	485
Length:	1400
Weight:	9
Length backrest:	500
Length seat:	500
Length legrest:	400
Weight capacity:	100

Subject to technical changes!



Date: 01.01.2024

Features:

- Perfectly suited for use in bathtubs or showers
- Sturdy cover made of "Sunsure" fabric
- Legrest and backrest continuously adjustable
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Guaranteed NO RUST
- Easy to clean
- 3 year warranty
- 5 year warranty on frame

Color selection Sunsure fabric

