

Datasheet: KSC 180 ASA FT PPK RP

Shower and Commode Chair for children

Price ex works Sargenroth:

Net list price [EUR]: 1.056,00

Technical data [mm/kg]:

Height:	850-950
Width:	520
Depth:	565
Weight:	9,5
Seat height:	420-520
Seat width:	400
Seat depth (without cushion):	415
Lower leg length:	300-400
Clearance under seat:	310-410
Drive-over width:	390
Weight capacity:	150

Subject to technical changes!



Date: 01.01.2024

Features:

- Height adjustment to 5 levels, allows to increase the seat height during the child's growing up
- Soft seat with toilet and hygiene opening
- Slide-out footrest
- Swing-out armrest (on the right, by request on the left)
- Backrest made of "Sunsure" fabric
- Removable backrest cushion with foam core and cover made of "Sunsure" fabric
- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm, all locking
- Pail with inside scale
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Padding thickness (unloaded) of backrest cushion: 50 mm
- Compression hardness of the foam: 3 kPa
- Easy to clean
- 3 year warranty
- 5 year warranty on frame

Color selection Sunsure fabric

