

Datasheet: KTLH 521

Hygiene and Isolation Cart

Price ex works Sargenroth:

Net list price [EUR]: 1.192,00

Technical data [mm/kg]:

Height:	1318
Width:	950
Depth:	505
Weight:	23
Height of top plate:	955
Height drawer size 1:	75
Width drawer size 1:	245
Depth drawer size 1:	355
Height drawer size 2:	175
Width drawer size 2:	245
Depth drawer size 2:	355

Subject to technical changes!



Date: 01.01.2024

Features:

- Liquid dispenser on slide rail (ESP), standard for half liter bottles (for liter bottles upon request)
- Drip tray for liquid dispenser
- Glove dispenser on slide rail (HSP)
- All plastic top plate
- 3 transparent drawers size 1
- 2 transparent drawers size 2
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Cart with Ø 75 mm twin-wheel swivel casters, 2 with lock
- Easy to handle docking system to attach hamper to cart
- Single hamper (capacity 120 liter/35 kg)
- Lid opening with foot pedal
- Hamper with Ø 75 mm twin-wheel swivel casters, 2 with lock
- Easy to clean
- 3 year warranty
- 5 year warranty on frame
- Delivery without pictured contents

Color selection lids

