

Datasheet: KTW 100 BEW

Blood sample cart

Price ex works Sargenroth:

Net list price [EUR]: 1.244,00

Technical data [mm/kg]:

Height:	1185
Width:	835
Depth:	390
Weight:	16
Height of top plate:	855
Height of drawer:	75 /175
Width of drawer:	245
Depth of drawer:	355

Subject to technical changes!



Date: 29.01.2024

Features:

- removable and padded armrest
- Holder for sharps disposable container
- waste container 10 liters
- Glove dispenser on slide rail (HSP)
- closed holder for skin antiseptic (max. 0.5 liter bottle)
- removable plastic shelf with one-sided upstand
- 4 transparent drawers, each with a load capacity of 8 kg
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm, 2 locking
- Easy to clean
- 3 year warranty
- 5 year warranty on frame
- Delivery without pictured contents