

# Datasheet: KTW 100 HS

## Hygiene and Isolation Cart

### Price ex works Sargenroth:

Net list price [EUR]: 622,00

### Technical data [mm/kg]:

Height:	1220
Width:	410
Depth:	385
Weight:	12
Height of top plate:	855
Height of drawer:	75 / 175
Width of drawer:	245
Depth of drawer:	355

Subject to technical changes!



Date: 01.01.2024

### Features:

- Liquid dispenser on slide rail (ESP), standard for half liter bottles (for liter bottles upon request)
- Glove dispenser on slide rail (HSP)
- All plastic top plate
- 4 transparent drawers, each with a load capacity of 8 kg
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm, 2 locking
- Easy to clean
- 3 year warranty
- 5 year warranty on frame
- Delivery without pictured contents