

# Datasheet: KTW 100 HS MRT

## MRI compatible Hygiene and Isolation Cart

### Price ex works Sargenroth:

Net list price [EUR]: 772,00

### Technical data [mm/kg]:

Height: 1220  
Width: 410  
Depth: 385  
Weight: 12,5  
Height of top plate: 855  
Height of drawer: 75 /175  
Width of drawer: 245  
Depth of drawer: 355

Subject to technical changes!



Date: 01.01.2024

### Features:

- Fulfills all safety requirements for stations up to 3 Tesla - MR CONDITIONAL
- Liquid dispenser on slide rail (ESP), standard for half liter bottles (for liter bottles upon request)
- Glove dispenser on slide rail (HSP)
- All plastic top plate
- 4 transparent drawers, each with a load capacity of 8 kg
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- 4 twin-wheel swivel casters (plastic/brass), Ø 75 mm, 2 locking, non-magnetic
- Easy to clean
- 3 year warranty
- 5 year warranty on frame
- Delivery without pictured contents
- For stations higher than 3 Tesla we offer customer specific product solutions, even all plastic equipment if required