

Datasheet: KTW 25 HS

Hygiene and Isolation Cart

Price ex works Sargenroth:

Net list price [EUR]: 622,00

Technical data [mm/kg]:

Height: 1220
Width: 410
Depth: 385
Weight: 12
Height of top plate: 855
Height of drawer: 175
Width of drawer: 245
Depth of drawer: 355

Subject to technical changes!



Date: 01.01.2024

Features:

- Liquid dispenser on slide rail (ESP), standard for half liter bottles (for liter bottles upon request)
- Glove dispenser on slide rail (HSP)
- All plastic top plate
- 3 transparent drawers, each with a load capacity of 8 kg
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm, 2 locking
- Easy to clean
- 3 year warranty
- 5 year warranty on frame
- Delivery without pictured contents