

# Datasheet: LTW 200

## Hygiene and Isolation Cart

### Price ex works Sargenroth:

Net list price [EUR]: 608,00

### Technical data [mm/kg]:

Height: 1455

Width: 705

Depth: 410

Weight: 19

Height of drawer: 165

Width of drawer: 530

Depth of drawer: 360

Subject to technical changes!



Date: 01.01.2024

### Features:

- 6 large, transparent drawers, each with a load capacity of 8 kg
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm, 2 locking
- Easy to clean
- 3 year warranty
- 5 year warranty on frame
- Delivery without pictured contents