

Datasheet: MC 231 S

Nursing Cart

Price ex works Sargenroth:

Net list price [EUR]: 1.090,00

Technical data [mm/kg]:

Height: 1425
Width: 1085
Depth: 505
Weight: 26
Usable shelf space: 410 x 330
Shelf spacing: 355
Capacity per shelf: 30
Subject to technical changes!



Date: 01.01.2024


Features:

- 3 removable shelves, size 1
- Shelves are removable without tools
- Roll-up front cover made of "Sunsure" fabric
- Liquid dispenser on slide rail (ESP), standard for half liter bottles (for liter bottles upon request)
- Glove dispenser on slide rail (HSP)
- Push-bar
- Single hamper
- Lid opening with foot pedal
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- 4 single-wheel swivel casters (plastic/stainless steel), Ø 100 mm, 2 locking
- Easy to clean
- 3 year warranty
- 5 year warranty on frame
- Delivery without pictured contents
- PLEASE NOTE: figure shows optional bags (WSW), not included in the price!

Color selection Sunsure fabric



Color selection bags



Color selection lids

