

**Net price ex works Sargenroth:**

Promotion price; no further discount [EUR]: 2.390,00

**Technical data [mm/kg]:**

Height:	1455
Width:	895
Depth:	525
Weight:	58
Height of working surface:	1060
Dimensions of working surface (WxD):	615 x 485
Interior height for modular baskets:	625
Subject to technical changes!	



Date: 18.01.2024

**Features:**

- Cabinet made of 19,6 mm RESOPAL® HPL plates with duroplastic film
- Stable frame made of RCN proprietary high impact plastic
- The material is scratch resistant, hygienic, UV stable, heat resistant and easy to clean
- Lockable drawer with 2 longitudinal dividers und 4 cross dividers
- Modular basket compartment (625 mm inner height) with 2 lockable swing doors
- 2 x ISO - modular basket 50, grey, 600x400x50 mm (including 1 longitudinal divider and 1 cross divider)
- 2 x ISO - modular basket 100, grey, 600x400x100 mm (including 1 longitudinal divider and 1 cross divider)
- 1 x ISO - modular basket 200, grey, 600x400x200 mm (including 1 longitudinal divider and 1 cross divider)
- Cabinet with 3-side upstand
- Push-bar on the left, height 1.110 mm (NOT PICTURED)
- Liquid dispenser on slide rail (ESP), standard for half liter bottles (for liter bottles upon request)
- Drip tray for liquid dispenser (not shown above)
- Glove dispenser on slide rail (HSP)
- Holder for sharps disposable container
- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 100 mm, 2 locking
- 3 year warranty
- 5 year warranty on frame

Color selection cabinet (MPC and MDC series)

