

Datasheet: MEW 50 Medical Cart

Price ex works Sargenroth:

Net list price [EUR]: 761,00

Technical data [mm/kg]:

Height: 875
Width: 425
Depth: 475
Weight: 17
Height of working surface: 875
Dimensions of working surface (WxD): 415 x 470
Interior height for modular baskets: 645
Subject to technical changes!



Date: 01.01.2024

Features:

- Frame made of RCN proprietary high impact plastic, sturdy and stable
- All plastic top plate
- Cover made of "Sunsure" fabric at the back and on both sides
- 2 x ISO - modular basket grey, 400x300x50 mm (including 1 longitudinal divider and 1 cross divider)
- 2 x ISO - modular basket grey, 400x300x100 mm (including 1 longitudinal divider and 1 cross divider)
- 1 x ISO - modular basket grey, 400x300x200 mm (including 1 longitudinal divider and 1 cross divider)
- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm, 2 locking
- Easy to clean
- 3 year warranty
- 5 year warranty on frame

Color selection Sunsure fabric

