

Datasheet: MPC 153 WSO

Multi-Purpose Cart

Price ex works Sargenroth:

Net list price [EUR]: 559,00

Technical data [mm/kg]:

Height: 1095

Width: 520

Depth: 440

Weight: 10,5

Usable shelf space: 410 x 330

Shelf spacing: 325

Capacity per shelf: 30

Subject to technical changes!



Date: 07.02.2024

Features:

- 3 removable shelves, size 1
- Shelves are removable without tools
- Push-bar
- Oval plastic wash bowl in red
- Washing bowl is removable (hung in the upper shelf)
- 7 liter capacity of the wash bowl
- Wash bowl dimensions: 450x300x130 mm
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- 4 single-wheel swivel casters (plastic/steel), Ø 100 mm, 2 locking
- Easy to clean
- 3 year warranty
- 5 year warranty on frame