

# Datasheet: MPC 153 WS MRT

## MRI compatible Multi-Purpose Cart

### Price ex works Sargenroth:

Net list price [EUR]: 688,00

### Technical data [mm/kg]:

Height: 1075  
Width: 520  
Depth: 440  
Weight: 10,5  
Usable shelf space: 410 x 330  
Shelf spacing: 325  
Capacity per shelf: 30  
Subject to technical changes!



Date: 01.01.2024

### Features:

- Fulfills all safety requirements for stations up to 3 Tesla - MR CONDITIONAL
- Plastic washbasin (integrated in shelf)
- 3 removable shelves, size 1
- Shelves are removable without tools
- Push-bar
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- 4 twin-wheel swivel casters (plastic/brass), Ø 100 mm, 2 locking, non-magnetic
- Easy to clean
- 3 year warranty
- 5 year warranty on frame
- For stations higher than 3 Tesla we offer customer specific product solutions, even all plastic equipment if required