

# Datasheet: MPC 353 ESP HSP

## Multi-Purpose Cart

### Price ex works Sargenroth:

Net list price [EUR]: 845,00

### Technical data [mm/kg]:

Height: 1360

Width: 930

Depth: 515

Weight: 16

Usable shelf space: 860 x 500

Shelf spacing: 315

Capacity per shelf: 30

Subject to technical changes!



Date: 01.01.2024

### Features:

- Liquid dispenser on slide rail (ESP), standard for half liter bottles (for liter bottles upon request)
- Glove dispenser on slide rail (HSP)
- 3 removable shelves, size 3
- Shelves are removable without tools
- Push-bar
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- 4 single-wheel swivel casters (plastic/steel), Ø 100 mm, 2 locking
- Easy to clean
- 3 year warranty
- 5 year warranty on frame